

Product List Serology 2010

| ABO System | | Clone | | | |
|------------------------------|---|--|--|--|--|
| Anti-A | | A11 H5 | | | |
| Anti-B | | B6 F9 | | | |
| Anti-AB | | A5 E10, B2 D7 | | | |
| Anti A1 | | | | | |
| Anti H (A2) | | 1093 4C 11 | | | |
| Rh System | | Clone | | | |
| Anti-D | IgM | MS 201 or RUM-1 | | | |
| Anti-D blend | IgM + IgG | TH28 / MS 26 | | | |
| Anti-D | IgG - incomplete | polyclonal | | | |
| Anti-C | | MS 24 or MS 273 | | | |
| Anti-c | | MS 33 or MS 35 | | | |
| Anti-Cw | incomplete | polyclonal | | | |
| Anti-CDE | IgM + IgG | MS24+MS201+MS26+MS80 | | | |
| Anti-E | | MS 80 / 258 or MS 12 / 260 | | | |
| Anti-e | | MS 16/21/63 or MS 62 / 69 | | | |
| Rh Control | Negative control for incomplete Rh sera | | | | |
| Kell System | | Clone | | | |
| Anti-Kell | | MS 56 or 601 | | | |
| Anti-Kell | coombsreactive | polyclonal, human IgG | | | |
| Anti-k (cellano) | coombsreactive | polyclonal, human IgG | | | |
| Anti-Kpa | coombsreactive | polyclonal, human IgG | | | |
| Anti-Kpb | coombsreactive | polyclonal, human IgG | | | |
| MNSs + P1 System | | Clone | | | |
| Anti-M | | M2-J4E-A7J or M11H2 | | | |
| Anti-N | | 1422 C 7 | | | |
| Anti-S | coombsreactive | polyclonal, human IgG | | | |
| Anti-s | coombsreactive | polyclonal, human IgG | | | |
| Anti-P1 | | 650 | | | |
| Duffy, Kidd and Lewis | | Clone | | | |
| Anti-Fya | coombsreactive | polyclonal, human IgG | | | |
| Anti-Fyb | coombsreactive | polyclonal, human IgG | | | |
| Anti-Jka | | MS 15 | | | |
| Anti-Jkb | | MS 8 | | | |
| Anti-Jka | coombsreactive | polyclonal, human IgG | | | |
| Anti-Jkb | coombsreactive | polyclonal, human IgG | | | |
| Anti-Lea | | LEA2 | | | |
| Anti-Leb | | LEB2 | | | |
| Rare Sera | | Clone | | | |
| Anti-Lua | coombsreactive | polyclonal, human IgG | | | |
| Anti-Lub | coombsreactive | polyclonal, human IgG | | | |
| Anti-Wra | coombsreactive | polyclonal, human IgM | | | |

Coombs Sera**Clone**

| AHG Serum | Coombs serum | Anti-IgG, Anti-C3d | | | |
|-----------|--------------|--------------------|--|--|--|

Supplements**Clone**

| | | | | | |
|--------------------|--------------------------------------|-----------------------|--|--|--|
| Bromelin | 200 testings | | | | |
| Bromelin Control | incomplete | polyclonal, human IgG | | | |
| Bovine Albumin 22% | | | | | |
| Bovine Albumin 30% | | | | | |
| Albu-Liss Reagent | | | | | |
| Alsever Solution | Cell stabilizer solution and diluent | | | | |